



## NEW: PAL System Vials & Vial Caps



- Borosilicate hydrolytic glass vials and polypropylene vials available
- Better sealing due to flat finish of contact surface to septum
- Improved bottom shape and thicker glass walls for fast and reliable transportation
- Reduced glass breakage through refined high quality glass
- Easy access to last drops due to slightly curved inner bottom shape
- Large labeling area for chemical resistant pencil inscriptions
- PAL System tested and approved

## Please see below all our available Vials and Vial Caps:

Vials	
PN CTC	Description
Vial-0.3-ND9-PP-100	Micro-Vial 2CV, 0.3mL Clear Polypropylene (PP), 12x32mm, with conical insert, fits ND9 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-1.5-ND9-AG-100	Vial 2CV, 1.5ml amber glass with Label. 12x32mm, 1st Class Hydrolytic Glass, fits ND9 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-1.5-ND9-CG-100	Vial 2CV, 1.5ml transparent clear glass with Label. 12x32mm, 1st Class Hydrolytic Glass, fits ND9 screw caps, Pk of 100 Pcs, screw caps not included
Vial-1.5-ND9-PP-100	Vial 2CV, 1.5mL Clear Polypropylene (PP), 12x32mm, with scale, fits ND9 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-10-DC20-CG-100	Vial 10CV, 10ml transparent clear glass with Label, for PAL Wash Modules. 46x22.5mm, DIN-Crimp Neck, Pk of 100pcs, w/o clip caps
Vial-10-ND18-AG-100	Vial 10CV, 10ml amber glass with Label. 46x22.5mm, fits ND18 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-10-ND18-CG-100	Vial 10CV, 10ml transparent clear glass with Label. 46x22.5mm, fits ND18 Screw Caps, PK of 100 Pcs, screw caps not included
Vial-20-ND18-AG-100	Vial 20CV, 20ml amber glass with Label. 75.5x22.5mm, fits ND18 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-20-ND18-CG-100	Vial 20CV, 20ml transparent clear glass with Label. 75.5x22.5mm, fits ND18 Screw Caps, Pk of 100 Pcs, screw caps not included

Caps	
PN CTC	Description
Cap-DC20-PE-SP13Y-100	Clip Cap for PAL Wash Modules. 22mm PE, Silicone/PTFE/Y-Slitted Septum 1.3mm, Pk of 100 pcs
Cap-ND18-St-SP15-100	Screw Cap 10CV and 20CV. ND18 Magnetic, Silver, Silicone white/PTFE blue Septum 1.5mm, Pk of 100 pcs
Cap-ND18-StC-SP20-100	Screw Re-Cap 10CV and 20CV. ND18 Magnetic, PTFE Metal Coating, Silicone/PTFE Septum 2.0mm, Pk of 100 pcs
Cap-ND9-St-SP10-100	Screw Cap 2CV. ND9, Magnetic, Gold, Silicone/PTFE Septum 1.0mm, Pk of 100 pcs
Cap-ND9-St-SP10S-100	Screw Cap 2CV. ND9, Magnetic, Gold, Silicone/PTFE/I-Slitted Septum 1.0mm, Pk of 100 pcs

Trial Set:	
Part No.	Description
PAL3-Vial-Sample-Set	Vial and Cap Sample Set consisting of: 2pc 1.5ml Clear and 2pc 1.5ml Amber Glass Vial, fits ND9 Screw Cap 2pc 10ml Clear and 2pc 10ml Amber Glass Vial, fits ND18 Screw Cap 2pc 20ml Clear and 2pc 20ml Amber Glass Vial, fits ND18 Screw Cap 4pc Screw Cap, ND9, Magnetic, Slitted Septa 1.0mm 4pc Screw Cap, ND18, Magnetic 4pc Screw Re-Cap, ND18, Magnetic

© 2017 CTC Analytics AG. All rights reserved. Neither this publication nor any part hereof may be copied, photocopied, reproduced, translated, distributed or reduced to electronic medium or machine readable form without the prior written permission from CTC Analytics AG, except as permitted under copyright laws. CTC Analytics AG acknowledges all trade names and trademarks used as the property of their respective owners.

PAL is a registered trademark of CTC Analytics AG | Switzerland

Date of print: 04.2017 | CTC Analytics AG, Industriestrasse 20, CH-4222 Zwingen, Switzerland, [www.palsystem.com](http://www.palsystem.com), Contact: [info@ctc.ch](mailto:info@ctc.ch)