

## LAB BENCH SERIES

Tables for heavy laboratory equipment like mass-specs and their accessories



## Instrument tables without integrated noise reduction cabinet



Sizes: 100x100x80, 150x100x80 and 200x100x80cm (WxDxH). Other sizes on request

Loading capacity: 500kg

Spacious cable duct with feedthroughs for cables and hoses

Mounted on high quality and lockable castor wheels

Height adjustable base for compensation of uneven floors

Feedthroughs for cables and tubing

Optional PC-cabinet available (see middle picture above)

## Instrument tables with integrated noise reduction cabinet



Fits all brands of mass-specs

Integrated noise reduction cabinet for one or two vacuum pumps

Noise reduction cabinet can be mounted to the left or to the right side of the table

Equipped with 30mm thick dampening material for maximum noise reduction

Pumps are protected by the Sonation APPS temperature and pump monitoring system

Integrated vibration decoupled oil tray on pull-out sliders for easy pump service

Sizes: 100x100x80, 135x100x80, 150x100x80 and 200x100x80cm (WxDxH). Other sizes on request

Loading capacity: 500kg

Spacious cable duct with feedthroughs for cables and hoses

Mounted on high quality and lockable castor wheels

Height adjustable base for compensation of uneven floors

Feedthroughs for cables and vacuumhose

## Electrically height adjustable side tables



Loading capacity: 150kg

Sizes: 60x80x765-1220, 120x80x765-1220cm (WxDxH). Other sizes on request

Integrated tray for waste bottles, samples etc.

Continuous height adjustment with 4 positions storable

Integrated pull-out drawer with soft stop

Cable duct with feedthroughs for cables, tubes etc.

Mounted multi-outlet power strip

Four lockable castor wheels with height adjustment for uneven floors (max. 20mm)

