


LAB TABLE SERIES HPLC

FLEXIBLE - MOBILE - ROBUST



Benefits

Flexible

- Many different widths, depths and heights available
- Electrically height adjustable
- Modular design

Mobile

- Ball bearing swivel castors with brakes
- Secure stand
- Very good running characteristics

Robust

- 250kg load capacity with even load distribution
- 150kg load capacity with uneven load distribution

Integrated cable duct

- Spacious design
- Velcro cable ties included

Integrated power distribution

- Socket strip integrated in the cable duct
- Multi-socket strips with different socket types available

Feedthroughs for cables and hoses

Fully compliant with laboratory standard 61010

Tested by an independent test laboratory

Heights

Fixed height models

60cm	(stocked)
70cm	(on request)
80cm	(stocked)
90cm	(on request)
100cm	(on request)

Height-adjustable models

53-83cm → 30cm travel	(stocked)
58-93cm → 35cm travel	(on request)
63-103cm → 40cm travel	(stocked)
68-113cm → 45cm travel	(on request)
73-123cm → 50cm travel	(on request)

Materials

Tabletop

- TRESPA® TOPLAB® PLUS
- Solvent resistant
- Easy to clean
- Excellent resistance to aggressive chemicals and cleaning agents
- Does not promote bacterial growth

Substructure and chassis

- Mechanically robust steel construction with powder-coated surface

table feet

- Aluminum



Widths and Depths

Widths

60cm	(stocked)
80cm	(on request)
100cm	(on request)
120cm	(stocked)

Depths

70cm	(on request)
80cm	(stocked)
90cm	(on request)
100cm	(on request)

Available accessories

Drawer module



Drawer module with soft-stop pull-outs.

Front panel height: 119mm
Drawer height: 68mm

Waste bottle tray



Tray with removable anti-slip mat, e.g. for waste bottles.

Depth: 23cm
Width: Table width - 20cm
Height: 5cm

Fixed floor plate



Shelf for e.g. a PC or accessories.
Material: TRESPA® TOPLAB® PLUS

IMPROVE YOUR SCIENCE

